




## CytoSoft®

Culturing Cells on an *In Vivo* Softness



# The Technology

Mechanotransduction refers to the processes through which cells sense and respond to mechanical stimuli by converting them to biochemical signals that elicit specific cellular responses. For example, changing to a stiffer substrate alters the differentiation potential of human mesenchymal stem cells to favor bone formation over cartilage and adipose tissues.

CytoSoft® plates have a thin layer of functionalized silicone with a stiffness between 0.2-64 kPa, allowing you to culture cells on a physiologically relevant matrix.


## Why CytoSoft®?

- Various areas throughout the body have significantly different degrees of matrix stiffness (ie. brain tissue ~0.2 kPa vs cartilage >64 kPa).
- Plastic and glass are >100,000X stiffer than physiological tissues.
- Culturing cells on an *in vivo* stiffness helps create *in vivo*-like results.
- Substrate stiffness affects propagation, gene expression, morphology, pluripotency, migration, differentiation, lineage specification




# Step 1: CytoSoft® Discovery Kit

## Discovering the Optimal Substrate Stiffness



The CytoSoft® Discovery Kit comes in a 7-pack of 6-well plates covering the entire physiological range of stiffness, including 0.2, 0.5, 2, 8, 16, 32 and 64 kPa.


- Seed cells across entire range in one experiment.
- Note differences in cells, as seen below.



As seen above, primary human dermal fibroblasts were cultured on 0.2, 8 and 64 kPa CytoSoft® plates. Fibroblasts grown on the 0.2 kPa and 64 kPa experienced high levels of F-actin stressfibers. Fibroblasts grown on the 8 kPa had an increase of adhesion sites and reduction of stressfibers.

## Step 2: CytoSoft® 6-Well Plates

### Further Research or Long-Term Cell Culture



- Buy an individual stiffness in a 5-pack of 6-well plates.
- Further identify cellular behaviors in a more narrow range of stiffness.
- Transition from long-term cell culture on plastic plates to CytoSoft® plates

---

The following cells have been cultured on CytoSoft® plates:

Adipose-Derived Stem Cell

Astrocyte

Fibroblast

HeLa

HUVEC

Cardiac fibroblast

Mesenchymal Stem Cell

MDCK

Myocyte

Neural Stem Cell

Glial cell

Cardiomyocyte

Longitudinal Smooth Muscle





## Step 3: CytoSoft® Imaging Plates

### High Resolution, Publication Quality Images



- #1.5 glass base
- Thinner silicone layer
- 24 and 96-well formats
- Obtain high resolution images



Live HeLa cells expressing fluorescent marker F-actin (lifeact-GFP). 60X 1.4NA oil objective. Scale bar 10 microns  
These images were taken on an 8 kpa *CytoSoft® Imaging* plate to demonstrate high resolution imaging capability

# Ditch the plastic plates and begin culturing cells on CytoSoft® today.

Product	Catalog #
CytoSoft® Discovery Kit	#5190-7EA
CytoSoft® 6-well Plate, Elastic Modulus 0.2 kPa	#5165-5EA
CytoSoft® 6-well Plate, Elastic Modulus 0.5 kPa	#5140-5EA
CytoSoft® 6-well Plate, Elastic Modulus 2 kPa	#5141-5EA
CytoSoft® 6-well Plate, Elastic Modulus 8 kPa	#5142-5EA
CytoSoft® 6-well Plate, Elastic Modulus 16 kPa	#5143-5EA
CytoSoft® 6-well Plate, Elastic Modulus 32 kPa	#5144-5EA
CytoSoft® 6-well Plate, Elastic Modulus 64 kPa	#5145-5EA
CytoSoft® Imaging 24-well Plate, Elastic Modulus 0.2 kPa	#5183-1EA
CytoSoft® Imaging 24-well Plate, Elastic Modulus 0.5 kPa	#5184-1EA
CytoSoft® Imaging 24-well Plate, Elastic Modulus 2 kPa	#5185-1EA
CytoSoft® Imaging 24-well Plate, Elastic Modulus 8 kPa	#5186-1EA
CytoSoft® Imaging 24-well Plate, Elastic Modulus 16 kPa	#5187-1EA
CytoSoft® Imaging 24-well Plate, Elastic Modulus 32 kPa	#5188-1EA
CytoSoft® Imaging 24-well Plate, Elastic Modulus 64 kPa	#5189-1EA
CytoSoft® Imaging 96-well Plate, Elastic Modulus 0.2 kPa	#5255-1EA
CytoSoft® Imaging 96-well Plate, Elastic Modulus 0.5 kPa	#5256-1EA
CytoSoft® Imaging 96-well Plate, Elastic Modulus 2 kPa	#5257-1EA
CytoSoft® Imaging 96-well Plate, Elastic Modulus 8 kPa	#5258-1EA
CytoSoft® Imaging 96-well Plate, Elastic Modulus 16 kPa	#5259-1EA
CytoSoft® Imaging 96-well Plate, Elastic Modulus 32 kPa	#5260-1EA
CytoSoft® Imaging 96-well Plate, Elastic Modulus 64 kPa	#5261-1EA


To order product or learn more, visit our website at

[www.AdvancedBioMatrix.com](http://www.AdvancedBioMatrix.com)

CytoSoft® plates require an extracellular matrix coating for cell attachment. Here are our recommended matrices:

PureCol® Type I Collagen #5005-100ML is the most widely used ECM for coating CytoSoft® plates.

- >99.9% Purity
- Cited in >2000 Publications
- Lot-to-Lot Consistency



---

Type III Human Collagen Solution, 10 mg	#5021-10MG
Type IV Human Collagen, Lyophilized, 5 mg	#5022-5MG
Fibronectin, Human, Lyophilized, 5 mg	#5080-5MG
Vitronectin, Human, Solution, 0.1 mg	#5051-0.1MG

# References

1. Modaresi, Saman, et al. "Deciphering the Role of Substrate Stiffness in Enhancing the Internalization Efficiency of Plasmid DNA in Stem Cells Using Lipid-Based Nanocarriers." *Nanoscale*, vol. 10, no. 19, 2018, pp. 8947–8952., doi:10.1039/c8nr01516c.
2. Wilson, Christina L. *In Vitro Models Of Brain For Study Of Molecular Mechanisms In Brain Disorder*. The University of Nebraska-Lincoln, 2016.
3. Wilson, C. L., Hayward, S. L., & Kidambi, S. (2016). Astroglial in a dish: Substrate stiffness induces astroglial in primary rat astrocytes. *RSC Advances*, 6(41), 34447-34457. doi:10.1039/c5ra25916a
4. Prager-Khoutorsky M, Lichtenstein A, Krishnan R, Rajendran K, Mayo A, Kam Z, Geiger B, Bershadsky AD. Fibroblast polarization is a matrix-rigidity-dependent process controlled by focal adhesion mechanosensing. *Nat. Cell Biol.* 2011; 13:1457–1465.
5. Gutierrez E, Tkachenko E, Besser A, Sundt P, Ley K, Danuser G, Ginsberg MH, Groisman A. High Refractive Index Silicone Gels for Simultaneous Total Internal Reflection Fluorescence and Traction Force Microscopy of Adherent Cells. *PLoS One*. 2011; 6:e23807.
6. Merkel R, Kirchgessner N, Cesa CM, Hoffmann B. Cell force microscopy on elastic layers of finite thickness. *Biophys. J.* 2007; 93:3314–23.
7. Schellenberg, A. et al. Matrix elasticity, replicative senescence and DNA methylation patterns of mesenchymal stem cells. *Biomaterials* 35, 6351–8 (2014).
8. Cesa, C. M. et al. Micropatterned silicone elastomer substrates for high resolution analysis of cellular force patterns. *Rev. Sci. Instrum.* 78, 034301 (2007).
9. Gutierrez, E. & Groisman, A. Measurements of Elastic Moduli of Silicone Gel Substrate with a Micro fluidic Device. *Plos One* 6 (2011).
10. Mori, H., Takahashi, A., Horimoto, A., and Hara, M. Migration of glial cells differentiated from neurosphere-forming neural stem/progenitor cells depends on the stiffness of the chemically cross-linked collagen gel substrate. *Neuroscience Letters*, Vol. 555, October (2013)
11. Banerjee, I., Carrion, K., Serrano, R., Dyo, J., Sasik, R., Lund, S. et al. Cyclic stretch of embryonic cardiomyocytes increases proliferation, growth, and expression while repressing Tgf- $\beta$  signaling. *J Mol Cell Cardiol.* 2015; 79: 133–144
12. Vertelov, G. et al. Rigidity of silicone substrates controls cell spreading and stem cell differentiation. *Sci. Rep.* 6, 33411; doi: 10.1038/srep33411 (2016).
13. Tkachenko E, Rawson R, La E, et al. Rigid Substrate Induces Esophageal Smooth Muscle Hypertrophy and EoE Fibrotic Gene Expression. *The Journal of allergy and clinical immunology*. 2016;137(4):1270-1272.e1. doi:10.1016/j.jaci.2015.09.020.
14. Sao, K. et al. Myosin II governs intracellular pressure and traction by distinct tropomyosin-dependent mechanisms. *Molecular Biology of the Cell* 30, 1170–1181 (2019).
15. Cooper, J. G. et al. Spinal Cord Injury Results in Chronic Mechanical Stiffening. *Journal of Neurotrauma* (2019). doi:10.1089/neu.2019.6540