



CERTIFICATE OF MATERIAL ORIGIN

Date: 10/22/2020

Revision: 03

Product: Collagen Membrane Film Discs, 23mm Diameter

| | |
|---|--|
| Supplier | Advanced BioMatrix, Inc. |
| Supplier Catalog No. | 5315-5EA |
| Supplier Location | California, USA |
| Name/Grade | Collagen Membrane Film, 23mm Diameter Discs |
| Animal Origin (including human) | |
| -Animal Species | Cephalopods |
| -Tissue Origin | Skin |
| -County of Origin | Pacific Ocean |
| -Country of Raw Material processing | Peru |
| Non-animal Source of Origin | |
| -Country of Origin | Not applicable |
| Animal Origin Material Adventitious Agent Risk Reduction Steps | Multiple washing and treatment steps |
| Method of Manufacturing | Extraction |
| Briefly describe Method of Manufacturing | Collagen membrane film is processed to render the material as a purified membrane from cephalopod skin via a multi-step process involving chemical soaking, extraction and pH treatment. |

The above information has been completed by an authorized representative of management and is complete and correct to my knowledge.

President

October 22, 2020

Corporate Headquarters
P.O. Box 502403
San Diego, CA 92150-2403

(800) 883-8220 (760) 929-0755
(510) 217-3452 Fax
www.advancedbiomatrix.com

Manufacturing Facility
5930 Sea Lion Place, Ste. 100
Carlsbad, CA 92010



Corporate Headquarters
P.O. Box 502403
San Diego, CA 92150-2403

(800) 883-8220 (760) 929-0755
(510) 217-3452 Fax
www.advancedbiomatrix.com

Manufacturing Facility
5930 Sea Lion Place, Ste. 100
Carlsbad, CA 92010