

# CERTIFICATE OF MATERIAL ORIGIN

Date: 03/26/25

Revision: 01

**Product:** Porcine, Type II Collagen, Solution

<b>Supplier</b>	Advanced BioMatrix, Inc.
<b>Supplier Catalog No.</b>	5465
<b>Supplier Location</b>	California, USA
<b>Name/Grade</b>	Porcine, Type II Collagen, Solution, ~3.0 mg/mL
<b>Animal Origin (including human)</b>	
-Animal Species	Porcine
-Tissue Origin	Cartilage
-County of Origin	United States
<b>Non-animal Source of Origin</b>	
-Country of Origin	Not applicable
<b>Animal Origin Material Adventitious Agent Risk Reduction Steps</b>	Multiple adventitious agent reduction and clearance steps
<b>Method of Manufacturing</b>	Purification
<b>Briefly describe Method of Manufacturing</b>	This product is highly purified preparation isolated from the triple helical domain of guanidine-HCl-extracted and pepsin-digested porcine joint cartilage collagen Type II.

The above information has been completed by an authorized representative of management and is complete and correct to my knowledge.



Senior Production Engineer

March 26, 2025